



5marfiecar


- *RF*
- *448KHz*
- *Upgraded*

Technical Features

MODEL	K-Smart Tecar
RF Frequency	300-448KHZ
Maximum Power	300W
Heads Size	20/40/50/60MM
Gross Weight	23kg
Package Dimension	50*48*42cm
Packaging	Alu Box+ Carbon Box

f Standard Accessories

Lastest Vesion of *Smart Tecar*



① RET Handle

② CET Handle

③ **4** CET Tips (20,40,50,60mm)

④ 3 RET Tips (20,40, 60mm)


⑤ Return pad

⑥ **Rotary switch** to adjust energy and time

⑦ Power switch

⑧ 10.4 inch touch screen

Real Pictures



Multifunctional Smart Tecar

01 In physiotherapy


Smart Tecar Therapy is a physical therapy that is part of the endogenous thermotherapy. It is able to significantly reduce recovery and rehabilitation time by stimulating the activation of all cellular repair processes. All this thanks to the application of radiofrequency that causes a thermal reaction in the tissues that stimulates the natural healing response with immediate anti-inflammatory and analgesic effects on muscles, tendons, cartilage or bones ligaments.

02 In beauty


Collagen requires a certain amount of heat to form. The heat produced by Smart Tecar in the tissue also stimulates collagen production. It helps to lift the skin and remove the wrinkle. At the same time, it also burns the fat and speeds up the elimination of body waste.

CET & I RET

The two systems, capacitive(CET) and resistive(RET), allow the operator to easily adjust the therapy according to the type of pathology.



Capacitive Technique: The energy emitted by Smart Tecar Concentrates in water-rich tissues: muscles, and organic fluids.



Resistive Technique: The energy emitted by Smart Tecar concentrates in the hard tissues: tendons, bone, and cartilage.



Advantages of Smart Tecar

01

In service

Factory direct supply
Professional and timely
aftersales service
Stable inventory and 72-hour
delivery time
Safe and punctual delivery
OEM & ODM service Multi-
color selection



Advantages of *Smart Tecar*

02 In machine

Double frequency

300KHZ and 448KHZ make RET and CET really different in deep and superficial penetration. The deeper penetration of RET can reach 10 CM without energy loss.

Multiple optional handles

Fascia set, D80mm handle, Foot pad, and T-shape handle.

Multi-language

Default: English and Spanish



Upgraded i

Upgraded i, lay the
recommende ; _c s for ten
body parts.

Unique, 二

Easier to a山侦『、h rotary
switch.


Larger screei

10.4 inch LED toi)reen

Higher power

300W output. Stable q transter.

Upgraded Interface



- Real-time parameters display
- One click to change the energy
 - ①Athermal (10%)
 - ②Medium thermal (30%)
 - ③Deep thermal (70%)



- Various handle options
- Each handle special for different body parts



- Ten body parts selection
- Each part with detailed recommended treatment protocols

Optional Accessories



01

Ø80mm

Do not need to work with return pad

Easier operation in

- > Elbow
- > Shoulder
- > Knee





**Do not need to work with
return pad**

**Easier operation in
> Cervical**

- > Spine
- > Palm
- > Foot



The heat flows through the foot to the whole body to speed up the blood flow and improve the metabolism.



Large size

- back

- thigh



Medium size

- calf
- shoulder




Small size

- acupuncture point

Functions

- > Physical immunity improvement
- > Muscular pain
- > Sport injuries
- > Myotenositis
- > Scar tissue
- > Sprains
- > Pelvic floor rehabilitation
- > Chronic pain
- > Facial management
- > Fat reduction
- > Postpartum repair



Usage (Guiding videos show the details)

Take "back treatment" as an example:

1. Input password " 888888 "
2. Apply RF cream to the b~ t
3. Put $t \angle \cdot \cdot \cdot$ re + ' "gainst the abdomen.
4. Choose the $\overline{\mathbb{E}}$ in the interface
5. Choose the treatme 、 p
6. Choose the treatment handle
7. Click on "Start"
8. Put the tip against the treatment area and move.

Notes

- * Need to work with RF cream, which helps energy transfer.
- * Need to work with return pad. (Except for the optional D80mm Hand and Foot pad)
- * Only use black CET tips carefully.
- * Keep the "is an 4 b i / En r U pad" (tip. The
- cIOSJ* is it s ir :e, " e tro : er he ie _ .
- * The silver RET tips need to be removed from the handle to clean after using.
- * Ask the feeling of the customer to avoid overheat because people have different feelings about the heat.



1. Do not put the return pad on the instrument.
2. Do not contact the working head with the return pad.
3. Do not touch the return pad lightly when the instrument is started.
4. The customer should completely touch the return pad to avoid electric shock.

Recommended Protocols

The treatment is divided into three periods.

- ① Athermal: 「博八刀」neration 5-10 mins.
 - ② Mediu. Thermal: Increase 沱 ie 07 F 乎.handle
slowly fr 2 , n 2
 - ③ Deep Thermal: Move the handle quickly to cover the larger treatment area with higher energy for 20 mins.
- The actual parameters should be set based on the conditions of the customer.

Treatment area	Treatment time	Interval	Sessions
Face (CET only)	30-35 mins	10 days	6 times
Abdominal fat reduction	45-60 mins	5-7 days	8 times
Finger	30 mins	Twice a week	6 times
Ankle	30-40 mins	Twice a week	8 times
Elbow	30-40 mins	Twice a week	6 times
Wrist	30-40 mins	Twice a week	6 times
Postpartum recovery	45 mins	Twice a week	12 times
Thigh	45 mins	Twice a week	8 times
Foot	40 mins	Twice a week	4 times
Caudal vertebra	40 mins	Twice a week	12 times
Shoulder	45-60 mins	Twice a week	8 times
Shoulder joint	35-40 mins	Twice a week	10 times
Arm	35-40 mins	Twice a week	8 times

Operation of standard handle ■ ■ CET

Picture



Treatment
area

Finger

Ankle

Elbow

Wrist

Tip

CET 20MM

CET 20MM

CET 40MM

CET 20MM

Operation

Athermal:

10% energy; fixed-point operation; 5-10 mins.

Medium Thermal:

20-30% energy; move the handle slowly; 10 mins.

Deep Thermal:

40% energy; move the handle quickly; 20 mins.

Interval: twice a week

Sessions: 6 times

Athermal: 10% energy; fixed-point operation for 5-10 mins.

Medium Thermal:

20% energy; move the handle slowly; 10 mins.

Deep Thermal: 30% energy; move the handle quickly; 20 mins. *Interval: twice a week*

Sessions: 8 times

Athermal: 10% energy; fixed-point operation for 5-10 mins.

Medium Thermal:

20% energy; move the handle slowly; 10 mins.

Deep Thermal: 30% energy; move the handle quickly; 20 mins. *Interval: twice a week*

Sessions: 6 times

Athermal: 10% energy; fixed-point operation for 5-10 mins.

Medium Thermal:

20% energy; move the handle slowly; 10 mins.

Deep Thermal: 30% energy; move the handle quickly; 20 mins. *Interval: twice a week*

Sessions: 4 times

Operation of standard handle -- RET

Picture



Treatment
area

Postpartu
m recovery

Thigh

Shoulder

Foot

Tip

RET 60MM

RET 60MM

RET 60MM

RET 40MM

Operation

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal:
20-30% energy; move the
handle slowly; 10 mins.

Deep Thermal:
40-60% energy; move the
handle quickly;
20 mins.

Interval: twice a week
Sessions: 12 times

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal: 20%
energy; move the handle
slowly; 10 mins. Deep

Thermal:
30% energy; move the
handle quickly;
20 mins.

Interval: twice a week
Sessions: 8 times

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal:
20-30% energy; move
the handle slowly;
10 mins.

Deep Thermal:
40-60% energy; move
the handle quickly;
20 mins.

Interval: twice a week
Sessions: 8 times

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal: 20%
energy; move the handle
slowly;
10 mins.





Deep Thermal: 30%
energy; move the handle
quickly;
20 mins.

Interval: twice a week
Sessions: 4 times

Operation of optional RET handle ■ ■ T-shape handle

Picture				
Treatment area	Caudal vertebra	Shoulder	Wrist	Foot
Tip	T-shape	T-shape	T-shape	T-shape
Operation	<p><u>Athermal:</u> 10% energy; fixed-point operation for 5-10 mins.</p> <p><u>Medium Thermal:</u> 20-30% energy; move the handle slowly; 10 mins.</p>			
	<p><u>Deep Thermal:</u> 40-60% energy; move the handle quickly; 20 mins.</p>			
	<p><u>Interval:</u> twice a week</p> <p><u>Sessions:</u> 12 times</p>			

Operation of optional RET handle ■ ■ 80MM handle

Picture				
Treatment area	Arm	Shoulder joint	Knee joint	Shoulder
Tip	80MM	80MM	80MM	80MM
Operation	<p><u>Athermal:</u> 10% energy; fixed-point operation for 5-10 mins.</p> <p><u>Medium Thermal:</u> 20% energy; move the handle slowly; 10 mins.</p> <p><u>Deep Thermal:</u> 30% energy; move the handle quickly; 20 mins</p> <p><u>Interval:</u> twice a week</p> <p><u>Sessions:</u> 8 times.</p>	<p><u>Athermal:</u> 10% energy; fixed-point operation for 5-10 mins.</p> <p><u>Medium Thermal:</u> 20-30% energy; move the handle slowly; 10 mins.</p> <p><u>Deep Thermal:</u> 40% energy; move the handle quickly; 20 mins.</p> <p><u>Interval:</u> twice a week</p> <p><u>Sessions:</u> 10 times</p>	<p><u>Athermal:</u> 10% energy; fixed-point operation for 5-10 mins.</p> <p><u>Medium Thermal:</u> 20-30% energy; move the handle slowly; 10 mins.</p> <p><u>Deep Thermal:</u> 40% energy; move the handle quickly; 20 mins.</p> <p><u>Interval:</u> twice a week</p> <p><u>Sessions:</u> 10 times</p>	<p><u>Athermal:</u> 10% energy; fixed-point operation for 5-10 mins.</p> <p><u>Medium Thermal:</u> 20-30% energy; move the handle slowly; 10 mins.</p> <p><u>Deep Thermal:</u> 40-60% energy; move the handle quickly; 20 mins.</p> <p><u>Interval:</u> twice a week</p> <p><u>Sessions:</u> 8 times</p>

Operation of optional RET handle-Fascia set

Picture



Treatment
area

Fasica set

Back

Shoulder

Calf

Acupuncture
point Small

Energy can only be 10%
in 5-10 seconds at a fixed
point. And then move
slowly for 7-10 mins.

Interval: twice a week
Sessions: 8 times

Large

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal:
20-30% energy; move the
handle slowly; 10 mins.

Deep Thermal:
40-60% energy; move the
handle quickly;
25 mins

Interval: twice a week
Sessions: 8 times.

Medium

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal:
20-30% energy; move
the handle slowly; 10
mins.

Deep Thermal:
40-60% energy; move
the handle quickly;
25 mins.

Interval: twice a week
Sessions: 8 times

Medium

Athermal: 10% energy;
fixed-point operation for
5-10 mins.

Medium Thermal:
20-30% energy; move
the handle slowly; 10
mins.

Deep Thermal:
40-60% energy; move
the handle quickly;
25 mins.

Interval: twice a week
Sessions: 6 times

Operation

Operation of optional CET handle--Foot Pad



Treatment area: foot

- 0 Athermal: 10-20% energy; 5-10 mins
- 0 Medium Thermal: 30-40% energy; 10 mins
- 0 Deep Thermal: 40-60% energy; 25 mins
- 0 Interval: 5-7 days
- 0 Sessions: 8 times

Precautions

- 0 Patients with impaired skin / neuropathy sensitivity.
- 0 Avoid treatment in patients with reduced skin sensibility in the area to be treated, unless the treatment is recommended by the patient's physician. If the treatment is to be used, test on the contralateral limb first to avoid tissue overheating.
- 0 Patients with a bacteriological inflammation
- 0 Diabetes mellitus: Patients with diabetes can often have neuropathy
- 0 Other conditions with reduced skin sensation
- 0 Patients coming in from the cold
- 0 Post-operative patients

Contraindications

- 0 Patients with pacemaker
- 0 Patients with bleeding gastrointestinal ulcers
- 0 Patients with an infusion pump and implants with electric lead
- 0 First 6 months after surgery - do not treat over the stomach and lower back (directly over the spine),
- 0 Treatment over the lower abdomen or for
- 0 Treatment of patients with a fever
- 0 Allergic reaction to cream / probe
- 0 Treatment of deep vein thrombosis
- 0 Treatment over the heart in patients with ischemic heart disease, local pulmonary embolism
- 0 Phlebitis
- 0 Bleeding areas where the skin has partial or open wounds

Kapha

Thanks

for your trust in

Kapha factory!



Jmapffecar